

Supplementary Document:**Detailed “Getting Started” Guide to 3D printing moulds for the SOL3D pipeline****1 1.0 BACKGROUND**

2 The field of 3D printing and bioengineering has developed rapidly during the last decades, in turn entering and
3 understanding a new field can be challenging. This document will provide more detailed insights to the
4 fundamentals of the techniques and their applications as well as limitations. The main document focusses on
5 the explanation of the developed pipeline, while this document illuminates the surrounding knowledge on 3D
6 printing and bioengineering.

7 Several microfabrication techniques have been adapted over the years to the needs of the biology research
8 community, and as a result the use of micropatterned substrates, microchamber devices and other engineered
9 substrates has increased exponentially. One of the most versatile combinations of techniques to obtain PDMS
10 devices for biological experiments is photolithography coupled with soft-lithography^{1,2}. Photolithography is
11 based on the deposition of layers of UV-sensitive photoresist of specified thickness, which is then exposed to a
12 UV source with either a photomask or by direct laser writing to create the desired design, before developing
13 the exposed photoresist. This process can be repeated for multiple layers and allows the creation of 2.5D
14 designs (i.e. multiple planar structures of different thickness stacked to form one single set of features)
15 (Diagram 1 A, B). While these techniques can be used to create advanced in vitro culture systems with micron-
16 scaled features, photolithography-based pipelines also present some limitations; for example, they can only
17 create a single layer at a time with a given height determined by the photoresist layer's properties and are
18 generally limited in the aspect ratio of the features that they can create. As a result, features are limited to
19 2.5D designs with defined thickness, lacking 3D volumes, curves, or interconnected shapes (Diagram 1B).
20 Moreover, generating multi-layer constructs for complex features involves multiple photolithographic steps,
21 which can be time-consuming, prone to errors, and costly (Diagram 1 A-B). While these limitations can be
22 obviated by recent improvements and optimisations in photolithography, such as grey scale photoresist^{3,4} and
23 high resolution 2-photon based lithographs^{5,6}, the result is often increased complexity in the fabrication
24 process, and in costs or availability of the necessary instruments. This creates the need for a technique that is
25 cost effective and provides features ranging from um to cm scale. Advancements in 3D printing, in price and
26 accessibility, render them as a great tool for biology labs. Here we outline the fundamental and commercially
27 available printer types, highlighting their advantages, disadvantages and estimated resolution. The printing
28 process overall always follow the same pattern though, with an idea, the design, the conversion to a 3D printer
29 suitable file format and the actual print. Later we go into detail, which software to use for the design and the
30 conversion to the correct file format, in this case tailored to the Phrozen 4k Mini.

31 Fused deposition modelling (FDM) allows to take a solid polymer, heat it to its melting point and extrude it
32 onto a printbed. The printer has a heat controllable hotend through which the filament is pressed, and a
33 coordinate system to move the hotend to a desired position in x, y and z. To construct an object, the design is
34 converted into g-code, giving the coordinates for the printer to move to, while extruding the polymer. This
35 type of printing takes rather long as every position is passed through and the object is constructed line by line.
36 However, it allows to use a vast range of different polymers with specific properties such as chemical stability
37 and durability. The resolution of this printer is determined by the filament size, which can range from 1,75mm
38 to 3mm, and the layer height, which can as low as 50- 100um but this ultimately depends on the printer
39 model.

40 Vat polymerisation is a specific type of 3D printing that generates constructs in a layer-by-layer fashion by
41 forming features onto a build plate using UV-curable resin. There are 3 main subtypes of vat polymerisation-
42 based 3D printers, defined by their UV light source, which dictate resin choice (due to light wavelength –

Supplementary Document:


Detailed “Getting Started” Guide to 3D printing moulds for the SOL3D pipeline

1 extensively reviewed⁷, as SLA printers usually use ~395nm light and LCD screens emit longer wavelength of
 2 ~405nm), part resolution, printing time, and printer cost. The original illumination technique for vat-
 3 polymerisation-based printing known as Stereolithography (SLA) utilised lasers, which function by scanning
 4 across designs on the build plate pixel-by-pixel to polymerise resin. These systems, such as the Formlabs 3B,
 5 offer sub-millimetre resolution (25µm X Y, 25 µm Z) and multiple proprietary optimised resin/curing options at
 6 the expense of printing speed, systems size, and cost, tending to be designed primarily for business use⁸.

7 The most common and cost-effective printer light sources are LCD screens, these illuminate each full layer of
 8 the design all at once, with X Y resolution dictated by the resolution of pixels within the screen, and Z
 9 resolution by the precision of the printer build plate motor. These printers offer similar resolution (35 µm X Y,
 10 µm Z) and, with the advent of mono-colour LCD screen technology, higher printing speeds than laser-based
 11 systems at a lower cost-of-entry and running⁹.

12 Digital Light Processing (DLP) printers employ a similar full-layer polymerisation method to LCD printers except
 13 designs are first projected onto a digital micromirror device (DMD) before being reflected onto the build plate.
 14 The benefit of DLP over LCD light sources is that the intensity of the illumination is uniform across the build
 15 area, meaning that µm-scale features are easier to manufacture due to a lack of pixel-pixel shadowing that can
 16 be present for LCD screens. However, DLP printers are more costly than LCD systems for equivalent resolution,
 17 and their resolution changes across print scales due to DMD manipulation of a single light source (51 µm X Y,
 18 10 µm Z) for USD 407 Anycubic Photon D2¹⁰.

19 Overall, driven by recent advancements of resolution and financial accessibility vat polymerisation 3D printers
 20 show the potential to bridge the gap between µm-resolution photolithography and mm-resolution fused
 21 deposition modelling (FDM) 3D printers (Diagram 1), whilst remaining economically accessible to any lab.
 22 Because of the positive characteristics UV vat polymerisation printers provide while remaining accessible,
 23 coupled with the rapid and easy fabrication of complex shapes in 3D, UV resin vat polymerisation represents in
 24 theory an ideal technique to create bespoke culture vessels, inserts and other devices to increase complexity
 25 within biological experiments without sacrificing control over culture conditions ^{11,12}.



Supplementary Document:

Detailed “Getting Started” Guide to 3D printing moulds for the SOL3D pipeline

1 While using PDMS as a material it is to note that, PDMS is a porous material that absorbs small molecules and
 2 growth factors, necessitating single usage of PDMS constructs and frequent preparation of new devices. For
 3 this reason, fast and reliable manufacturing of these devices is often a crucial limit step in several experimental
 4 pipelines.





A**B****C****D**

Diagram 1: Description of 3D vat polymerization and soft lithography and scale comparison across methodologies

(A) Schematic overview of the UV resin vat polymerisation process and post-processing with a time estimate for production. Printing times range depend on the printed volume but for cell culture device range from 10min to 1h. (B) Schematic overview of photolithography and an estimated time to completion of a print. Print time is an estimate for feature creation across the whole silicon wafer and alters in between prints, depending on the design. Overall, UV resin vat polymerisation is faster compared to photolithography, as it uses whole field illumination and no single point illumination as photolithography. (C) Scale comparison of 3D printing methods UV resin vat polymerisation, FDM and photolithography. Comparison of maximal resolution achieved with UV resin vat polymerisation, with current printers, and photolithography. (D) Representation of feature designs achievable with UV resin vat polymerisation and photolithography. UV resin vat polymerisation offers a wide range of feature dimensions and complexity compared to photolithography which underlies technical limitations for designs.

Supplementary Document:

Detailed “Getting Started” Guide to 3D printing moulds for the SOL3D pipeline

1 2.0 DESIGN

2 When designing a part, a clear idea of the final construct and its function are important to make the process
 3 quick and easy. For example, visualising the spatial relationships between a 3D printed mould and the PDMS
 4 cast it is intended to form in the printing – can you print the features you need? Casting – does the resin you
 5 plan to print with enable PDMS curing with/without enamel coating? Demoulding – does the part have
 6 multiple/complex components? Can it be demoulded in a single step? and Cell culture function – is it a stencil
 7 device? Does it need to be thin for imaging? all dictate the design strategy at the CAD level.

8 Software

9 There are several CAD software’s online, both subscription based and open source, the most common are
 10 [Fusion360](#) (subscription based) and [Tinkercad](#) (open source). Although both are offered by Autodesk Inc.
 11 Tinkercad is a high-level design software, useful for creating designs rapidly from a menu of shapes. Whilst
 12 Fusion360 is a more powerful software which builds designs from fundamental sketch components, enabling
 13 the creation of more complex and adaptable designs. Although Fusion360 is subscription based, for those in
 14 research institutions with a yearly renewable education license can be obtained so it is beneficial to opt for
 15 Fusion360 over tinkercad in most cases. In addition, there is a wealth of design tutorials on YouTube to help
 16 get to grips with Fusion360.

17 Picking a resin and printer (SLA)

18 Depending on the intended function of the printed part, resin and printer selection is vital to the success of the
 19 construct. For example, ABS like resin has good mechanical properties but poor compatibility with PDMS
 20 casting. Prints with small features <200um and curved surfaces should be printed with high-resolution resin to
 21 improve fine feature formation and curved surface smoothness. Additionally, it is important to consider the
 22 size and number of parts to be fabricated as printers have different sizes and shapes of build plate. A
 23 breakdown of the capabilities of the resins used in the Serio lab and printers they are suitable for is below.

Resin	Layer thickness settings	Type of prints	Currently printed on
4K aqua Gray	25um, 50um	High res, curved	4K, Mars 3
Elegoo ABS like	50um (from online repository)	Strong prints, non PDMS cast	4K, Mars 3, Photon S
Next Dent Ortho Clear	50um (default)	Biocompatible/non leaching – non PDMS cast	4K
Flexible X	50um (non-optimised)	Flexible prints/ non PDMS casting	4K
Premium tough	50um (non-optimised)	Parts with good mechanical properties, non PDMS cast	4K

Supplementary Document:

Detailed “Getting Started” Guide to 3D printing moulds for the SOL3D pipeline

1

2

3 Designing parts for PDMS casting

4 As mentioned, the design of parts intended for PDMS casting requires a clear plan of the final function of the
 5 part. Here, it is important to understand the limits of printing resolution, PDMS mechanical properties, and
 6 practical aspects of the casting process such as using a positive mould to fabricate a negative with PDMS.

7 For example, whilst parts with features of 100 µm in XYZ and spacing can be manufactured, moulds for
 8 microgrooves of 100um width, depth and spacing require tolerances of 50um to be added to the spacing of
 9 grooves. This accounts for bleaching that occurs between the channels as resin drainage is worse for these
 10 types of features than for standalone features. Less spacing than this between channels results in shallower
 11 features than specified at this depth. Another factor to note is that when printing micro features (<200um)
 12 decreasing layer thickness does not guarantee improved resolution. For squared/single profile features, having
 13 less layers by using 50um layer thickness can give better results than with 25um layer thickness for the same
 14 print due to potential overcuring of the smaller layers. Lower layer thicknesses are helpful for prints with
 15 features that have more complex shapes e.g., round/v shaped microgrooves as smaller layers smooth
 16 gradients, providing a more accurate replication of the CAD with impacts on cell behaviour when plated in
 17 constructs.

18 Minimum feature sizes optimised:

19 4K resin on 4K/Mars 3 printer @25um layer

Shape	CAD Dimension (um)	Actual Dimension (um) XYZ
3 Dimension (feature)	200x200x200, 200 spacing	230x220x240, 150 spacing
2 Dimension (grooves)	100x250, 200 spacing	120x >100x 240, 190 spacing

20

21 Additionally, both printing and demoulding tolerances should be accounted for when designing a part for
 22 PDMS casting. As a rule, an outer wall thickness of 2mm is the minimum suggested to enable successful
 23 printing of larger (taller >5mm) moulds and to prevent mould destruction during cast demoulding. For internal
 24 design spacing, at >1mm between any features and the walls is recommended (more if you can to help
 25 prevent damaging casts when demoulding), with >500um between features >500um in height. Aspect ratios
 26 are important here so take this into consideration when designing too. See diagram below.

27


28

29

30

Supplementary Document:

Detailed “Getting Started” Guide to 3D printing moulds for the SOL3D pipeline




1

2 Designing parts in Fusion360

3 In fusion 360, all designs are generated from parametrically constrained **sketch** objects. This ensures that
4 designs can be easily modified by altering any parameters in the sketch, without having to restart the design
5 process. Sketches are formed into 3D objects through **extrusion**. Modification of 3D object can then be
6 performed with several operations including: **filet**, **chamfer**, and **shell** to generate a final part. Below is a short
7 example of how to construct a stencil mould design as demonstrated in the SOL3D manuscript.

8 Stencil Device Tutorial

9 Select **Sketch** from the ‘solid’ menu, click on a plane, and create the base of the stencil device mould by
10 clicking again to create the sketch, typing in the desired length of each dimension. Here, lengths of 40mm X
11 and 30mm Y were chosen. Once happy, click ‘Finish Sketch’ on the right-hand tool bar to complete the sketch.




12

13 Once the sketch is complete you can **extrude** the shape into a 3D base object by right clicking on the
14 completed sketch and selecting **extrude** from the menu. The Sketch will become darker, and another definable
15 dimension icon will appear on the sketch. Here the sketch was extruded 1mm.


Supplementary Document:

Detailed “Getting Started” Guide to 3D printing moulds for the SOL3D pipeline



1

- After completing the base section, repeat the **sketch** and **extrude** functions to create the 3D structure that will become the mould for the PDMS casts (below). All sketches so far were made with the ‘centre rectangle’ tool that automatically snaps to the centre of the sketch plane, simplifying the design process. Keep in mind the final dimensions of the cast you want to make here as they will be based on the size of the hole in this extrusion NOT the extrusion itself.




7

- The dimensions chosen here (above) were 27mm X and 7mm Y and 4mm Z to account for losses when using the **shell** tool on the upper face of the extrusion. 1mm shell (below) leaves a 1mm wall on each side, removing 2mm from both X and Y dimensions, leaving a hole 25 mm X 5 mm Y and 4 mm Z.

Supplementary Document:

Detailed “Getting Started” Guide to 3D printing moulds for the SOL3D pipeline




1

- Now the basic mould is complete you can customize the wells you would like to print. In this example we will design a mould to cast a stencil with 4 equally sized 3.5 mm X 1 mm Y 2mm Z wells, with 3.5 mm X 3 mm Y 2 mm Z funnels above, spaced 2 mm apart, with a tolerance of 2.5 mm on each side for easy demoulding. The benefit of SLA printing here is that you can customize well size, shape, and configuration to experimental needs, and editing the Fusion360 design file is quick and easy through rolling back the **history marker (red)**.
- To begin the design of the wells, start from the top (the funnel). First, create a sketch as before of the top side of the funnel (3.5 mm X, 3 mm Y). Since the base was generated using a ‘centre rectangle’ sketch the base centreline forms an easy reference point to align sketches too (the tool will also snap to this point to assist you). Here I also sketched a **construction line** (2.5mm from the wall), selected from the ‘sketch palette’ on the right of the screen, which does not form part of the final design but can help in aligning sketches on a design.


Supplementary Document:

Detailed "Getting Started" Guide to 3D printing moulds for the SOL3D pipeline



1

- 2 After creating the funnel sketch you can **extrude** as before to the 2mm Z length specified in our design
3 requirements. To form the tapered funnel shape, use the **chamfer** tool by selecting the two long edges
4 (3.5mm) at the top of the extrusion and inputting 1 mm as the distance. This will constrain the well size to 1
5 mm in Y, the desired dimension. Next, set the angle of the **chamfer** by selecting ‘distance and angle’ in the
6 option menu. This will vary depending on the depth of the well, and the size of the funnel and well. In this case
7 an angle of 62° was used to maximise the funnel size with continuous angle (Below). Fusion360 will not let you
8 input an angle that is too high as it over-constrains the object.




9

10

Supplementary Document:

Detailed “Getting Started” Guide to 3D printing moulds for the SOL3D pipeline


- 1 After finishing the **chamfer** for the funnel, simply **extrude** the well 2 mm to match the height of the mould walls (Below). All there is to complete is copying the well within the mould, and copying the whole mould across the base.



- 4
- 5 To copy parts of your designs, use the **rectangular pattern** tool (below). In the options panel on the right be sure to select ‘faces’, before selecting all the faces you want to pattern. Here you can see there are 9 faces selected, this is higher than you might expect because you need to select the small faces at the bottom of the funnel **chamfer**. This is a caveat of the **chamfer** tool where you are unable to create a flush surface as doing this over-constrains the face. Once all the faces are selected you change from ‘objects’ to ‘directions’ in the side menu and select on the design an edge (selected in blue below) that is parallel to the desired pattern direction. From here you set the number of wells (4) and the spacing between them (5.5 mm) to give a 2 mm gap between adjacent wells and a 2.5 mm gap to the mould wall. Spacing is set from the starting point of the selected objects in the desired direction, which in this case includes the 3.5 mm X dimension of the well. This is important to take into consideration when designing different well sizes.


Supplementary Document:

Detailed “Getting Started” Guide to 3D printing moulds for the SOL3D pipeline



1

- 2 See below the completed single cast mould after using the **rectangular pattern** tool to copy the well within the mould.




4

- 5 To complete the design, use the **rectangular pattern** tool again to copy the mould you just completed across the full area of the base. The difference between this patterning and the previous is that the mould is copied both above and below its original position, requiring the selection of **Symmetric** from the ‘direction type’ option in the menu panel on the right (below).


Supplementary Document:

Detailed “Getting Started” Guide to 3D printing moulds for the SOL3D pipeline



1

2




3

- 4 Above is the completed PDMS stencil device mould, complete with 12 3.5 mm X 1 mm wells spaced 2 mm apart, and a 2 mm deep funnel on top to facilitate easy cell seeding – concluding the design tutorial.
- 5 Once a design is complete it needs to be exported to an stl. (Stereolithography mesh) to be sliced and read by the 3D printer. This can be done easily in Fusion360 with File>export>.stl (below).

Supplementary Document:

Detailed “Getting Started” Guide to 3D printing moulds for the SOL3D pipeline



1

- 2 NB: Be sure to continually save your designs as you go, Fusion360 does not autosave. Conveniently, however,
3 each save point can be recovered so if certain alterations don't work you can always go back. Rolling the
4 history marker back also helps to modify designs without starting from the beginning.
- 5 If you have any doubts about how to use Fusion360 don't forget there are several online sources of tutorials
6 such as the Product design online [channel](#). These tutorials can also enable you to design more complex, multi-
7 part constructs to suit your individual needs.
- 8

9 3.0 SLICING

10 Softwares


- 11 Once a design is finalised and you're ready to print you need to process the stl. Design file into slices to be
12 printed layer-by-layer by the printer. This is achieved with a slicing software. There are several different slicers
13 available, however the most used are Chitubox for SLA printers and Cura for FDM printers.

14 Chitubox

- 15 Chitubox is an easy-to-use slicing software with support for most SLA printers and resins built in, and a host of
16 basic model manipulation functions to help get the best print possible. After importing a design, Chitubox
17 allows the manipulation of orientation, scale, positioning, cloning of parts. As well as some more invasive
18 functions such as shelling parts or drilling holes in designs. See below the window set up of Chitubox 1.9.1.

Supplementary Document:

Detailed “Getting Started” Guide to 3D printing moulds for the SOL3D pipeline



1

2 To correctly set a print going using Chitubox there are 3 things to consider.

3

- Print orientation/supporting
- Balancing the build plate
- Resin settings

6 Setting up a printer for the first time

7 If you have not used a printer before you need to set it up in the chitubox ‘settings’ tab. As chitubox is the main supplier of boards for most SLA printers they support a wide range of printers and resins with manufacturer published default resin settings. You need to select ‘add a printer’ and scroll to the system you have. For the Serio lab this will either be the Phrozen sonic mini 4K, or the Elegoo Mars 3. You can have any printers you use saved in your ‘settings’ and add/remove more as you please.

12 NB: New printers and resin settings are updated with chitubox updates, make sure the software is as up to date as possible.


14 NBB: As chitubox software updates roll through sometimes printers stop supporting the sliced file type. If this happens it will most likely be fixed by a printer firmware update. These are slightly tedious but will be explained in detail on manufacturer websites if required.

17 Print orientation/supporting

18 Once you have added your printer you and imported a part by dragging in your .stl file you need to consider the orientation of your part, and whether it needs supporting. To do this you can use the ‘move’, ‘rotate’, ‘scale’, and ‘mirror’ tools to change the location, scale, or orientation of your print – see below. Print orientation becomes extremely important for printing parts with microfeatures such as grooves where features are acutely impacted by factors such as pixel shape and orientation, grooves are printed at 90°.


Supplementary Document:

Detailed “Getting Started” Guide to 3D printing moulds for the SOL3D pipeline



1

- Once you are happy with the size and position of your print, you need to consider whether you need to support it or not. As prints are constructed layer-by-layer it is important to identify and support parts of your designs that overhang or are not flat as these areas cannot be printed without additional supports. When supporting parts, it is also good practice to angle designs to limit the risk of failure of overhangs, although this will reduce print resolution in some cases – lower layer thicknesses will help negate these losses, at a trade off with print time. See below a print requiring angle changes and supports.




8

- Here it is important to note that the LED array has rectangular pixels, so isotropic designs print better when aligned with the array. Meaning that regular elongated pattern such as grooves should be printed at 90°, compared to round, dot-like features which appear the best at 0°.


Supplementary Document:

Detailed “Getting Started” Guide to 3D printing moulds for the SOL3D pipeline

A



B



Supplementary Document:

Detailed "Getting Started" Guide to 3D printing moulds for the SOL3D pipeline

1 **Balancing the build plate**

2 After finalising the supports for a print, you need to check that the build plate is balanced. This is important as
3 unbalanced prints will result in differential pulling force on the FEP film and build plate when retracting from
4 each polymerised layer, increasing the risk of FEP rupture but also compromising print resolution. For
5 large/single prints it is important to centre the part on the build plate with the move tool. For multiple prints it
6 is good practice to organise the build plate in a symmetrical fashion, with the tallest parts in the centre.

7 **Resin settings**

8 Finally, you need to select the appropriate resin settings for your chosen printer, resin, and layer thickness
9 from default or optimised parameters in the 'settings' tab. Understanding these parameters can help in fine
10 tuning print quality for a specific print, or for optimising a new printer/resin.


11 In short, the settings can be grouped into 3 types: exposure, light off delay, and lifting speed/distance. With
12 'bottom exposure time' dictating how well adhered prints will be to the build plate. 'Exposure time' dictating
13 the layer thickness (how much resin is polymerised for each specified 'layer height'). And the coupling of lifting
14 speed/distance with 'light off delay' dictating how far and fast the build plate is moved up and down between
15 each layer, and how long the LCD waits before re-illuminating. In the lab, our experience has found that
16 'exposure time' has the most significant impact on print quality variation, and that slower lifting speeds (and
17 associated light off delay) will increase the quality of prints (especially for higher viscosity resins) at the expense
18 of print speed. 'Transition layers' are a series of layers with decreasing exposure time between the bottom
19 layers and the normal layers, providing better adhesion of the top and bottom layers of your print. Whilst they
20 can help limit the 'elephant foot' phenomena that occurs from wider resin polymerisation at the base of prints
21 from the high exposure, they are not compulsory. A repository of the current optimised print parameters for
22 resins can be found at the end of this section. A light off delay calculator can also be found in the github
23 repository to help if modifying settings from default.

24 NB: New updates to chitubox can sometimes alter the structure and parameters in resin settings so it is
25 important to keep track of changes in case they alter PDMS compatibility or print quality.


26

Supplementary Document:

Detailed “Getting Started” Guide to 3D printing moulds for the SOL3D pipeline



- 1 All you need to do now is hit the ‘slice’ button and save your sliced .ctb file to a virus free USB stick. Once in the slicing screen you can also see an estimate of the consumable’s usage of each print, the print time, and a reconstruction of each printed layer. This reconstruction is a good opportunity to double check for any unsupported islands on prints. There is also a brief summary of the resin settings at the bottom of the window which is also helpful to check you have selected the correct parameters.



- 7
8
9

Supplementary Document:

Detailed “Getting Started” Guide to 3D printing moulds for the SOL3D pipeline

1 4.0 PRINTING

2 Preparing to print

3 Cleaning

4 Before starting a print, it is important to check the printer is clean, this will help ensure print quality, reduce
5 the chance of failures, and increase the longevity of the system. Regardless of printer type there are several
6 important components to be checked before starting to print:

7 Build plate

8 Check there is no residual resin on the build plate from previous prints and clean off. If left unchecked resin
9 can be partially cured from exposure to natural light when the printer is open and can result in poor print
10 adhesion to the build plate. If particularly worried, use IPA to help remove resin.

11 Resin VAT and LCD Screen

12 VAT – **EMPTY** – if the printer hasn’t been used in a while, or you are printing in a different resin, ensure you
13 clean the inside and outside of the resin VAT thoroughly. Do this by wiping down the exterior (with or without
14 IPA) and gently wiping the top and bottom of the FEP film with a tissue (softer/finer than paper towel, also
15 with or without IPA) to remove any dust/residual resin.

16 VAT – **FULL** – If the printer is being used with the resin already present in the VAT it is important to check there
17 is no solid debris left in the resin from previous prints. If there are, this will likely indent or rupture the FEP film
18 during the next print, potentially causing damage to the printer from resin leaking from the VAT into
19 components. Check for resin debris by gently running the mixing spatula through the resin (this will also mix
20 the resin): feeling for increased resistance – parts polymerised to the FEP film, or for solid objects that appear
21 at the surface when moving through the resin. If debris are found ensure to remove them and check
22 thoroughly that there are no more present. For parts on the FEP film, try to remove them with gentle pressure
23 using the spatula (be careful not to damage the FEP) and repeat the check for debris. The best way to remove
24 all debris is to change the resin in the bath.

25 Changing the resin

26 To change SLA resins, you need to drain and clean the resin VAT completely before adding the new resin. If
27 there is a lot of resin left in the bath (more than 1cm depth) it is best to drain this into a ‘dirty’ resin bottle so it
28 can be reused (do not mix different resins – have a separate bottle for each). If there is only a little resin left in
29 the vat it is likely contaminated and not worth keeping. To clean waste resin all printers, have a vat cleaning
30 function. This function uniformly illuminates the full build area for a user defined period, polymerising a layer
31 of resin to allow it to be peeled off the FEP film with ease and can be disposed of safely – leaving a clean vat
32 behind. A normal exposure is 20 seconds and will allow polymerisation without fusion of the resin to the FEP
33 film. A caveat of the vat cleaning function is that the screen is by design smaller than the area of the FEP film,
34 so additional mopping of resin in the surrounding area is required. A tip here to prevent marking the FEP is to
35 conduct the vat cleaning but mop excess resin from around the now-polymerised resin before peeling off the
36 semi-polymerised resin. Once the vat is clear of resin you can add your desired resin to the bath.

37

Supplementary Document:

Detailed “Getting Started” Guide to 3D printing moulds for the SOL3D pipeline

1

2

3 Mixing the resin

4 Before printing all resins need to be **mixed** properly to ensure even distribution of photoinitiators and other
5 resin components such as colour pigments. Without proper mixing, print quality will be varied and risk of
6 failure is higher.

- 7 • Adding resin from the bottle:
8 ○ The best way to mix resin to ensure good mixing whilst avoiding creating bubbles is to gently
9 invert the resin bottle 5-10 times before pouring into the vat.
10 • Mixing resin in the vat:
11 ○ Use a plastic spatula to gently mix the resin without touching the FEP film. Look out for
12 changes in resin colour and distribution of pigments as an indicator of mixing.

13 Getting rid of bubbles

14 In many cases, in particular with more viscous resins, bubble formation during mixing is unavoidable. Bubbles
15 in the resin can result in print failure so it is very important to remove as many of them as you can before
16 printing. The best way to remove bubbles from the resin is to let the resin sit for 10-15 minutes before starting
17 a print to allow them to surface and disappear. If you are short on time the tip of the spatula (or razor/ pipette
18 tip) can be used to gently remove bubbles at the surface of the resin manually.

19 Securing the resin bath and build plate

20 The final step before setting a print going is to check that the resin vat is firmly secured to the printer by
21 checking the screw on each side is tight. It is important for print resolution to have the vat as tightly secured as
22 possible to the printer to ensure even distribution of forces over the FEP film on retract and uniform exposure
23 of resin across the bottom of the bath.

24 Setting a print going

25 Once you have the printer cleaned and resin mixed, turn the printer on, plug in the USB, find the file you would
26 like to print, and press go!

27

28

29 5.0 POST PROCESSING

30 Post processing of resin prints serves 2 purposes, to **wash off** excess resin left after printing, and to **finish the**
31 **UV polymerization** of the resin in its final design.


32 Removing prints from the build plate

33 After print completion and before post processing can begin you need to remove the part from the build plate.
34 Depending on the resin you printed with there will be excess present on the build plate, remove this with the

Supplementary Document:

Detailed “Getting Started” Guide to 3D printing moulds for the SOL3D pipeline

1 plastic spatula back into the resin vat before removing the build plate from the print. Once most of the excess
 2 resin is removed you can detach the print from the build plate using a combination of a razor blade and the
 3 plastic spatula – see below:



4

5

Prize up corner
with razor

6

Push razor further
under print


Slip spatula underneath
and push rest of print off

7 It is important to prize up the part **gently**, firstly to prevent scratching or indenting the build plate which might
 8 impact the quality of subsequent prints, but also to limit the risk of destroying the not-fully-polymerised part.

9 See below the optimized protocol for all resin prints.

10 **Post processing – all resins:**

- 11 1. 10min in Sonicator in IPA (dirty wash) – in instrument room with -80 freezer
- 12 2. 10min in wash of wash & cure in IPA (clean wash)
- 13 3. Cure 1h (wash and cure system)



14

15 10 mins sonicator (dirty wash)

10 mins wash (clean wash)

16 1 hour UV curing – ensure parts
are completely dry of IPA before
curing

- 17 - For biocompatible resin use **fresh** isopropanol not filtered
- 18 - If printing very small <250 µm features, less/no sonication can be performed to preserve the delicate
19 features. If extremely worried, conduct manual washing with IPA, although this will impact the curing
20 of PDMS if casting is required

Supplementary Document:

Detailed “Getting Started” Guide to 3D printing moulds for the SOL3D pipeline

- 1 - **IMPORTANT – DO NOT** leave prints in IPA for prolonged periods as overexposure can result in warping. Warping is worst from the sonicator as the IPA here gets hot (>50°C after use).
- 2
- 3



Supplementary Document:

Detailed "Getting Started" Guide to 3D printing moulds for the SOL3D pipeline

1 Deciding which pipeline to use

We developed two strategies for casting PDMS in 3D printed moulds, which you can find highlighted in the main manuscript Figure 1G. The main difference between these two protocols is the additional coating with enamel paint after the post curing process. The two main contributing factors to decide which pipeline to use, are the size of the features and the type of mould. As a rule of thumb, complex and closed moulds with little air contact, such as sandwich moulds used for figures 8 and 7 will require the additional coating. Additionally, if your designs do not contain small scale features <100um, it is worth it to immediately use the paint protocol shortening the troubleshooting process. But this decision has to be taken for each design individually. We added here some design examples that take these points into consideration

10

Coated pipeline	Non-coated pipeline
 <p>Muscle device mould: A complex mixture of small and large-scale features in a sandwich mould. There is limited air contact, a large contact area and complex features that make the demoulding itself challenging. The mould has been enamel coated, even in the small features, through angle variation while air brushing.</p>	 <p>Microwell mould: This mould has a large surface area but only small features. Here, we used the non-coating pipeline, so that the features remain at their original size.</p>
 <p>Plating device mould: In this case, we decided to use the painting protocol, as no feature will be covered and the thick layer of PDMS can cure better. This might be also possible with the non-coating protocol</p>	 <p>Groove mould: This mould has a large surface area but only small features. Here, we used the non-</p>


Supplementary Document:**Detailed “Getting Started” Guide to 3D printing moulds for the SOL3D pipeline**

but for better demoulding we decided to use the coating pipeline, so that the features remain the original size.	coating pipeline, so that the features remain the original size.
--	--

1

2

- 3 To equally distribute the enamel paint on the surface airbrush system is of great help and if it is carefully done
- 4 and the print rotated multiple times from different angles, one can reach even into deep and complex features.
- 5 In our experience if paint does not reach the area, the design is too complex and the PDMS will not demould
- 6 properly even if it would cure.



Supplementary Document:**Detailed “Getting Started” Guide to 3D printing moulds for the SOL3D pipeline**

- 1
- 2
- 3 1. Soft lithography for micro- and nanoscale patterning.
- 4 <https://www.nature.com/articles/nprot.2009.234.pdf>.
- 5 2. Qin, D., Xia, Y. & Whitesides, G. M. Soft lithography for micro- and nanoscale patterning.
- 6 *Nat Protoc* **5**, 491–502 (2010).
- 7 3. Rammohan, A. *et al.* One-step maskless grayscale lithography for the fabrication of 3-
- 8 dimensional structures in SU-8. *Sens. Actuators B: Chem.* **153**, 125–134 (2011).
- 9 4. Waits, C. M., Modafe, A. & Ghodssi, R. Investigation of gray-scale technology for large
- 10 area 3D silicon MEMS structures. *J. Micromechanics Microengineering* **13**, 170 (2003).
- 11 5. Think big. Print nano. Your partner for high-precision additive manufacturing.
- 12 <https://www.nanoscribe.com/en/>.
- 13 6. UpNano. <https://www.upnano.at/technology/#scale-applications>.
- 14 7. Gong, H., Beauchamp, M., Perry, S., Woolley, A. T. & Nordin, G. P. Optical approach to
- 15 resin formulation for 3D printed microfluidics. *RSC Adv.* **5**, 106621–106632 (2015).
- 16 8. High Resolution SLA and SLS 3D Printers for Professionals | Formlabs.
- 17 https://formlabs.com/uk/?&utm_source=google&utm_medium=cpc&utm_campaign=UK-EMEA-Prospecting-Search_Brand-Trademark-Brand-EN-Exact-Paid-Adwords&utm_term=formlabs%20printer&utm_content=formlabs_printers&utm_device=c&_bt=345203181553&_bk=formlabs%20printer&_bm=e&_bn=g&_bg=70222759115&gclid=Cj0KCQjwsp6pBhCfARIsAD3GZuauRLTauENV3-KUHR17vNNe0U026sjAqTjInAnjqoYeaubTq4E4TEcaArxJEALw_wcB.
- 23 9. Phrozen Sonic Mini 4K Resin 3D Printer | Phrozen Technology: Resin 3D Printer
- 24 Manufacturer. <https://phrozen3d.com/products/sonic-mini-4k-resin-3d-printer-phrozen>.
- 25 10. Anycubic Photon D2 – ANYCUBIC-UK. <https://uk.anycubic.com/products/anycubic-photon-d2?variant=44183762272541>.
- 27 11. Balakrishnan, H. K. *et al.* 3D Printing: An Alternative Microfabrication Approach with
- 28 Unprecedented Opportunities in Design. *Anal Chem* **93**, 350–366 (2021).
- 29 12. Kamei, K. *et al.* 3D printing of soft lithography mold for rapid production of
- 30 polydimethylsiloxane-based microfluidic devices for cell stimulation with concentration
- 31 gradients. *Biomed. Microdevices* **17**, 36 (2015).
- 32 -